

IWC
PARIS
2008



5th International Whey Conference



UNIVERSITÄT
HOHENHEIM



WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION

Fermentation and enzymology as tool for functional whey ingredients/products

Lutz Fischer

**Institute of Food Science and Biotechnology
Department of BIOTECHNOLOGY
University of Hohenheim, STUTTGART, Germany**



WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION

Sweet cheese whey (components in g L⁻¹)

Water
(938-943)

pH 5.9-6.6

Total Solids
(57-62)

Carbohydrates (47)
-Lactose (47)

Milk Fat (0.5)
- Triglycerides (0.25)
- Phospholipids (0.15)
- Diglycerids (0.05)
- Free Fatty Acids (0.05)

Proteins (6.2)
- *b*-Lactoglobulin (3.0)
- *a*-Lactalbumin (1.2)
- Immunoglobulin-G (0.7)
- Proteose Peptone (0.6)
- Serum Albumin (0.4)
- Other Protein (0.3)
 # Immunoglobulin-A (0.05)
 # Lactoferrin (0.045)
 # Lactoperoxidase (0.025)
 # Lysozyme (0.002)

Minerals (4.5)
- Potassium (1.5)
- Chloride (1.1)
- Phosphorous (0.7)
- Calcium (0.6)
- Sodium (0.5)
- Magnesium (0.1)
- Trace Elements (0.03)
 # Zinc (0.0015)
 # Iron (0.0006)
 # Iodine (0.0005)
 # Copper (0.0002)

NPN (0.132)
- Urea (0.08)
- Free Amino Acids (0.025)
- Cholin (0.015)
- Orotic Acid (0.012)

Vitamins (0.075)
- Vitamin B5 (0.004)
- Vitamin B2 (0.0015)
- Vitamin C (0.0015)
- Vitamin B6 (0.0005)



Fermentation using whey




WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION


Comparison: typical
chemically defined medium

component	amount/L
glucose	10 g
NH ₄ Cl	1 g
K ₂ HPO ₄	0.5 g
MgSO ₄ • 7 H ₂ O	0.2 g
FeSO ₄ • 7 H ₂ O	0.01 g
CaCl ₂ • 2 H ₂ O	0.01 g
NaCl	0.5 g
trace element solution	1 ml

„molecular formula“ of a bacterial cell :



Whey is an undefined
natural complex medium



The quality and/or quantity
of the **carbon/energy**
source and the
nitrogen source is an issue


- *b*-galactosidase activity
- additional nitrogen source(s)



WP as fermentation medium



Processed whey concentrates (deproteinized)



sophorose lipid



Single cell protein (SCP)



- ~ yeasts (GRAS)
*K. marxianus**
- $Y_{X/S}$ 0.48-0.52
- ~ feed
- ~ (food)

Single cell oil (SCO)

- ~ yeasts, fungi
Cryptococcus curvatus
- ~ triglycerides (27 g L⁻¹)
 - è N-limitation
 - è unsaturated FA
(ca. 50% oleic acid)
 - è possible use: biodiesel



Biosurfactants (SL)

- ~ yeasts
Candida boidinii
- ~ biotensid (420 g L⁻¹)
 - è on cell debris from
C. curvatus



Constituent ^a	Final concn (mM) ^a	
	MPL _R ^c	MPL (amt/liter)
Lactose	58.4	58.4 (20 g)
Sodium thioglycolate*	4.38	4.38 (0.5 g)
Tween 80* (% vol/vol)	0.1	0.1 (1 g)
K ₂ HPO ₄ · 3H ₂ O*	7	
Na ₂ HPO ₄		2.11 (0.3 g)
NaH ₂ PO ₄ · H ₂ O		5.8 (0.8 g)
Potassium acetate*	15	75 (7.35 g)
Ammonium citrate*	10	
Orotic acid*	3.2	
Pyruvic acid*	0.9	
Formic acid*	0.2	
Oleic acid*	0.033	
NaCl	10	
CaCl ₂ *	5	2.5 (0.27 g)
MgSO ₄ · 7H ₂ O*	0.8 ^d	0.8 (0.2 g)
MnSO ₄ · H ₂ O	0.12 ^d	0.12 (0.02 g)
FeSO ₄ · 7H ₂ O (μM)	70 ^d	0.23 (0.066 mg)
Micronutrients* (μM)		
(NH ₄) ₆ Mo ₇ O ₂₄ · 4H ₂ O	0.003 ^d	0.003 (0.0037 mg)
CoCl ₂ · 6H ₂ O	0.03 ^d	0.03 (0.007 mg)
H ₃ Bo ₃	0.4 ^d	0.4 (0.025 mg)
CuSO ₄	0.01 ^d	0.01 (0.0016 mg)
ZnSO ₄ · 7H ₂ O	0.01 ^d	0.01 (0.0029 mg)
Bases*		
Adenine	0.37 ^d	0.37 (0.05 g)
Guanine	0.33 ^d	0.33 (0.05 g)
Xanthine	0.32	0.32 (0.05 g)
Uracil	0.45 ^d	0.45 (0.05 g)
2'-Deoxyguanosine*	0.19 ^d	
Hypoxanthine*	0.37	
Amino acids		
L-Alanine	2.24 ^d	1.12 ^d (0.1 g)
L-Arginine	0.29 ^d	1.84 ^d (0.32 g)
L-Asparagine	1.9 ^d	2.27 (0.3 g)
L-Aspartate	7.5	3.75 ^d (0.5 g)
L-Cysteine · HCl	2.28 ^d	1.6 (0.28 g)
Glycine	2.66 ^d	2.13 ^d (0.16 g)
L-Glutamate	5.9 ^d	2.25 ^d (0.38 g)
L-Glutamine	2 ^d	2 (0.3 g)
L-Histidine	3.22 ^d	1.3 ^d (0.2 g)
L-Isoleucine	3.8 ^d	2.75 ^d (0.36 g)
L-Leucine	1.52 ^d	4.57 ^d (0.6 g)
L-Lysine · HCl	0.27 ^d	2.4 ^d (0.44 g)
L-Methionine	3.35 ^d	0.8 ^d (0.12 g)
L-Phenylalanine	1.21 ^d	2.06 ^d (0.34 g)
L-Proline	4.7	8 ^d (0.92 g)
L-Serine	1.9 ^d	3.42 ^d (0.36 g)
L-Threonine	0.42 ^d	2.52 ^d (0.3 g)
L-Tryptophan	1 ^d	0.5 ^d (0.1 g)
L-Tyrosine	1.1 ^d	0.66 ^d (0.12 g)
L-Valine	4.27 ^d	4 ^d (0.48 g)
L-Cystine	0.42 ^d	
Vitamins (μM)		
B ₁₂ (cobalamin)*	7.4 ^d	7.4 (0.01 mg)
Riboflavin*	2.66 ^d	0.53 (0.2 mg)
Calcium pantothenate*	2 ^d	0.2 (0.1 mg)
Niacin*	14 ^d	1.4 (0.17 mg)
Folic acid*	0.72 ^d	0.072 (0.032 mg)
Pyridoxal*	2.45 ^d	0.245 (0.05 mg)
Pyridoxamine*	0.83	
p-Aminobenzoate*	2.92 ^d	
Biotin*	0.04 ^d	
Thiamine*	3 ^d	
L-Ascorbic acid* (mM)	2.8	
Spermine · 6H ₂ O*	2.5	

- ~ lactic acid bacteria need a huge number of supplements (auxotrophic)
- ~ typical complex medium: **MRS** (peptone, meat extract, yeast extract, glucose, salts, minerals)
- ~ **production of *L. casei*** using supplemented and diluted, partly deproteinated whey^a (from powder)
 - è supplements: ammonium salts (2.5 g L⁻¹), yeast extract (0.25 g L⁻¹)
 - è biomass yield to lactose: $Y_{X/S} = 0.165 \text{ kg kg}^{-1}$



by-products: **lactic acid** and **EPS**



Lactic acid production in whey



q production of pure L(+)- or D(-)-lactic acid*

- 80.000 t anno⁻¹, 90% by LAB fermentation
- polylactic acid: biodegradable polymers
- synthon in chemistry
- *food*: preservative, flavour agent, acidity adjuster

Lactic acid isomers produced by common dairy lactobacilli

Name of bacteria	D(-) Lactic acid	L(+) Lactic acid	DL (±) mixture
<i>L. acidophilus</i>	No	No	Yes
<i>L. delbrueckii</i> subsp. <i>lactis</i>	Yes	No	No
<i>L. delbrueckii</i> subsp. <i>bulgaricus</i>	Yes	No	No
<i>L. helveticus</i>	No	No	Yes
<i>L. casei</i>	No	Yes	No
<i>L. paracasei</i> subsp. <i>tolerans</i>	No	Yes	No
<i>L. paracasei</i> subsp. <i>paracasei</i>	No	Yes ^a	No
<i>L. rhamnosus</i>	No	Yes	No
<i>L. plantarum</i>	No	No	Yes
<i>Lc. lactis</i>	No	Yes	No
<i>S. thermophilus</i>	No	Yes	No
<i>Leuconostoc</i> sp.	Yes	No	No

+ yeast extract

pH 5.5, 37°C,
4 g L⁻¹ h⁻¹

continuous production using immobilised cells (entrapment by gels)

--> **productivity**
12 g L⁻¹ h⁻¹

^a Some strains produce D(-) and L(+) lactic acid (Source: Curry and Crow, 2003a).



Catalysis on lactose

WHEYVOLUTION

WHEYVOLUTION

WHEYVOLUTION

WHEYVOLUTION

WHEYVOLUTION

Lactose-free milk, whey products

- Increased digestibility, sweetness, solubility



enzymatic hydrolysis

Lactose (Gal-Glu)

trans galactosylation

Galactosyl oligosaccharides (GOS)

- Increased solubility
- Prebiotic (baby foods)
- Laxative



Vivinal® GOS

enzymatic catalysis ?

Lactulose (Gal-Fru)

- Most popular laxative
- 10-times value of lactose
- Prebiotic (FOSHU list)
- Baby foods as bifidogenic factor
- Increased sweetness, solubility





Synthesis using a hydrolase?




WHEYVOLUTION

WHEYVOLUTION

WHEYVOLUTION

WHEYVOLUTION

WHEYVOLUTION



q reversed hydrolysis
(thermodynamically controlled)

§ a_w -value

§ concentration donor/acceptor



q transgalactosylation
(kinetically controlled)

§ nucleophilic character

§ concentration of acceptors



geothermal heated
marine sediments
(Stetter, Uni Regensburg)



Quelle: ArchiMeDes

- q *Pyrococcus furiosus* (extremophilic species of Archaea)
- q enzyme CelB:
 b-glucosidase as well as *b*-galactosidase activity
- q stability: $t_{1/2}$ (100°C) = 85h, $t_{1/2}$ (110°C) = 13h
- q CelB produced in *E. coli* (1 mkat L⁻¹ = 16.6 Mill. U L⁻¹)


Unusual biocatalyst



Expanded bed reactor system



WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION



- 1: column, optimum temperature
- 2: immobilised enzyme (Cel B)
- 3: substrate solution (processed whey)
- 4: water bath
- 5: pump
- 6: product solution

$$\text{Expanded Bed} \Rightarrow \frac{H_{\text{total}}}{H_{\text{beads}}}$$

- covalently immobilised Cel B
- whey UF permeate
- at 75°C



Continuous lactulose synthesis



WHEYVOLUTION

WHEYVOL

WHEYVOL

WHEYVOL

WHEYVOL

Whey-UF-permeate	sugar	H_{total}/H_{beads}			
		start t_0	1.3	1.6	2.0
diluted + 0.1 M Frc	Lactose [g/L]	45.61	5.93	8.20	11.02
	hydrolysis [%]		87.0	82.0	75.8
	Lactulose [g/L]		0.96	1.25	1.44
	yield [%]		2.4	3.3	4.2
+ 0.5 M Frc	Lactose [g/L]	44.17	6.68	8.49	10.95
	hydrolysis [%]		84.9	80.8	75.2
	Lactulose [g/L]		5.66	6.47	6.83
	yield [%]		15.1	18.1	20.6
concentr. + 0.1 M Frc	Lactose [g/L]	109.66	34.11	41.56	48.05
	hydrolysis [%]		68.9	62.1	56.1
	Lactulose [g/L]		2.36	1.71	1.75
	yield [%]		3.1	2.5	2.8
+ 0.5 M Frc	Lactose [g/L]	105.13	20.81	33.65	39.61
	hydrolysis [%]		80.2	68.0	62.3
	Lactulose [g/L]		3.84	10.62	10.62
	yield [%]		4.6	14.9	16.2



WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION

Selective protein hydrolysis



- ~ whole cell fermentation
- ~ commercially available enzymes (pancreatic proteases, etc.)
- ~ proteases from dairy starter cultures

Swiss, Finland:

Evolus




Japan: **Calpis**



dairy products with
ACE inhibitors



L. helveticus: proteolytic enzymes



Native-PAGE of the soluble proteolytic enzymes from *Lactobacillus helveticus*: active-staining (H-Ala-Pro-pNA; H-Ala-pNA) and silver-staining

(A: active-staining, B: silver-staining; Phast™-Gradient-Gel 8-25, GE Healthcare Bio-Science, Uppsala, Sweden).



Substrate specificity



WHEYVOLUTION

WHEYVOLUTION

WHEYVOLUTION


WHEYVOLUTION

WHEYVOLUTION


	relative activity [%]	
	whole cell	soluble enzymes
H-Ala-Pro-pNA	100*	100**
H-Ala-pNA	8.9	67
H-Pro-pNA	2.9	3.4
H-Val-pNA	0.67	3.2

*100% = 81 nkat/mL


**100% = 60 nkat/mL




Captopril



Lisinopril



Isoleucine-Proline-Proline (IPP)



Valine-Proline-Proline (VPP)

Examples of drug ACE-inhibitors
[Estler and Schmidt, 2007]

Structural formula of the bioactive, anti-hypertensive tripeptides produced in skim milk



Sequences of whey proteins

WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION
WHEYVOLUTION

a-Lactalbumin B (123 amino acids)

EQLTKCEVFRELKDLKGYGGVSLPEWVCTTFH
TSGYDTEAIVENNQSTDYGLFQINNKIWCKNDQ
DPHSSNICNISC DKFLNNDLTNNIMCVKKILDK
VGINYWLAHKALCSEKLDQWLCEKL

b-Lactoglobulin B (162 amino acids)

LIVTQTMKGLDIQKVAGTWYSLAMAASDISLLD
AQSAPLRVYVEELKPTPEGDLEILLQKWENGEC
AQKKIIAEKTKIPAVFKIDALNENKVLVLDTDYK
KYLLFCMENSAPEQSLACQCLVRTPEVDDDEAL
EKFDKALKALPMHIRLSFNTPQLEEQCHI



- ~ whey and its various formulations (permeates, powders, concentrates) can be regarded as versatile raw materials for further refinement by microbial cells or/and selective enzymes
- ~ there is and probably will not be „the one big process“ but a wide variety of bright biocatalytic applications to produce new functional valuables
 - è food industry
 - è cosmetic industry
 - è pharmaceutical industry
 - è chemical industry
- ~ entrepreneurial decisiveness is demanded from the dairy industry in order to introduce new methods and expertise and to expand their interdisciplinary activities



Karin Jaindl, Timo Stressler, Thomas Eisele,
Anja Manns, Sabine Lutz-Wahl, Bertolt Kranz,
Wolfgang Claaßen

Prof. Dr.-Ing. Jörg Hinrichs and his group

... and the financial support by the BMBF
- AiF-FV 14787 (via FEI, Bonn)
- BioChance PLUS-3 FKZ 0313137